


# Magnum<sup>®2</sup> Knotless Fixation Implant

Minimal size. Maximum bone fixation.



The OPUS<sup>®</sup> Magnum<sup>2</sup> combines strong bone fixation and an increased bone lock surface area that maximizes pull out strength.

- Knotless design allows for speed, strength and a more repeatable repair
- The butterfly wing design provides superior bone lock purchase and a strong sub-cortical fixation
- Helps prevent the possibility of implant "pull-back" by setting bone lock prior to suture lock
- Suture is locked within the implant reducing the possibility of suture slip, gap formation or repair failure
- TensionLock™ allows for lateralization of the cuff to restore the anatomic footprint with fine .5mm adjustments



Magnum<sup>2</sup> Knotless  
Implant

OPUS AutoCuff® System  
Perfecting the Art of  
Rotator Cuff Repair



SmartStitch® with PerfectPasser®  
Fast, Easy, Reliable



Magnum<sup>2</sup>  
Knotless Fixation Implant




MagnumWire®  
High Tensile Strength Suture



Completed Rotator Cuff Repair  
with OPUS Magnum<sup>2</sup>

## Magnum<sup>2</sup> Implant Specifications



Part Number	Description
OM-1502CS	Magnum <sup>2</sup> Implant, Case
OM-9010	Magnum <sup>2</sup> Drill
OM-9015	Magnum <sup>2</sup> Drill Guide and PathFinder® Set
OM-9020	Magnum <sup>2</sup> Punch
OM-9120S	Magnum <sup>2</sup> Procedure Set PerfectPasser Polyester (Mixed) Suture Cartridges
OM-9121	Magnum <sup>2</sup> Procedure Set PerfectPasser MagnumWire (Mixed) Suture Cartridges
OM-9126	Magnum <sup>2</sup> Plus PerfectPasser MagnumWire (White) Suture Cartridge
OM-9127	Magnum <sup>2</sup> Plus PerfectPasser MagnumWire (Co-Braid) Suture Cartridge
OM-9029S	Magnum <sup>2</sup> Procedure Set M-Connector Polyester (Mixed) Suture Cartridges
OM-9021	Magnum <sup>2</sup> Procedure Set M-Connector MagnumWire (Mixed) Suture Cartridges
OM-9026	Magnum <sup>2</sup> Plus M-Connector MagnumWire (White) Suture Cartridge
OM-9027	Magnum <sup>2</sup> Plus M-Connector MagnumWire (Co-Braid) Suture Cartridge

To order or for more information, please contact your ArthroCare Sports Medicine representative today or call 800-797-6520.



**ArthroCare Sports Medicine**  
680 Vaqueros Avenue  
Sunnyvale, CA 94085-3523  
main phone 408-736-0224  
order entry **800-797-6520**  
order entry fax 888-994-2782  
[arthrocaresportsmedicine.com](http://arthrocaresportsmedicine.com)

**ArthroCare Europe AB**  
Skeppsbron 2  
111 30 Stockholm  
Sweden  
phone +46 8 546 172 00  
fax +46 8 546 172 39  
[info@arthrocare.se](mailto:info@arthrocare.se)

**CAUTION:** Federal (USA) law restricts this device to sale by or on the order of a physician.  
© 2009 ArthroCare Corporation. All rights reserved. ArthroCare, OPUS, AutoCuff, Magnum, PerfectPasser, SmartStitch, Pathfinder, and MagnumWire are trademarks and/or registered trademarks of ArthroCare Corporation.