

# OPUS® AutoCuff® Magnum® X Knotless Fixation Implant with Independent Tensioning


Be Independently Inspired



Magnum X with Independent Tensioning combined with the SmartStitch® Suture Passing System provides a transosseous tunnel equivalent for an ideal repair in an efficient manner.

- Magnum X allows the surgeon to tension each suture limb individually, or together, to independently apply and maintain even pressure on the suture limbs and repair construct for ideal contact
- Knotless design locks securely under the cortex, reducing the possibility of gap formation
- TensionLock™ allows for lateralization of the cuff to restore the anatomic footprint with fine .5mm adjustments
- New snare shuttles simplify suture loading
- Achieves a bone lock with unparalleled pullout strength

Magnum X provides its greatest advantage in the lateral row of a double row configuration. If the suture limbs are different lengths, Magnum X allows the surgeon to independently tension them and provide equally distributed pressure over the restored footprint.




Magnum Knotless  
Suture Bridge

## The AutoCuff System

The Ultimate In  
Rotator Cuff Repair



SmartStitch with PerfectPasser™  
Fast, Easy, Reliable




SmartStitch instantly places  
an Incline Mattress® Stitch



MagnumWire®  
The Strongest Suture




Magnum M  
Medial Implant for Superior  
Double Row Repairs



OPUS Knotless Double  
Row Repair

## Magnum X Implant Specifications



| Part Number | Description  |
|-------------|--|
| OM-1503CS   | Magnum X with Independent Tensioning, Case   |
| OM-9033     | AutoCuff Magnum X Procedure Set with MagnumWire  |
| OM-9036     | AutoCuff Magnum X Plus Knotless Implant Set with MagnumWire (White) Suture Cartridge                 |
| OM-9037     | AutoCuff Magnum X Plus Knotless Implant Set with MagnumWire (Co-braid) Suture Cartridge              |
| OM-9131     | AutoCuff Magnum X Procedure Set with MagnumWire  |
| OM-9136     | AutoCuff Magnum X Plus Implant Set with PerfectPasser MagnumWire (White) Suture Cartridge            |
| OM-9137     | AutoCuff Magnum X Plus Implant Set with PerfectPasser MagnumWire (Co-braid) Suture Cartridge         |
| OM-1525CS   | Magnum M Fixation Implant, Case  |
| OM-9531     | Double Row Procedure Set with Magnum M, Magnum X, PerfectPasser and 4-Runner™ Suture Shuttle (Cases) |
| OM-1000CS   | 4-Runner Suture Shuttle, Case  |
| OM-9010     | Magnum Drill (3.2mm)   |
| OM-9015     | AutoCuff Drill Guide and Pathfinder® Set   |

To order or for more information, please contact your ArthroCare Sports Medicine representative today or call 800-797-6520.



**ArthroCare Sports Medicine**  
680 Vaqueros Avenue  
Sunnyvale, CA 94085-3523  
main phone 408-736-0224  
order entry **800-797-6520**  
order entry fax 888-994-2782  
**arthrocaresportsmedicine.com**

**ArthroCare Europe AB**  
Baggensgatan 25  
111 31 Stockholm  
Sweden  
phone +46 8 546 172 00  
fax +46 8 546 172 39  
**info@arthrocare.se**

**CAUTION:** Federal (USA) law restricts this device to sale by or on the order of a physician.

© 2009 ArthroCare Corporation. All rights reserved. All trademarks, service marks and logos listed herein are registered and unregistered trademarks of ArthroCare Corporation, its affiliates or subsidiaries or a third party who has licensed its trademarks to ArthroCare.